Camera Module CM 900®

Zubehör zu HAAG-STREIT Spaltlampen
Accessoire pour lampes à fente HAAG-STREIT
Accessory to HAAG-STREIT slit lamps

Gebrauchsanweisung
Mode d'emploi
Instruction manual

Vorwort

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Introduction

We would like to thank you for your decision to purchase this HAAG-STREIT product.

If the instructions in this manual are carefully followed we are confident that this product will give you reliable and trouble-free usage.


Purpose of use

The **Camera Module CM 900®** is an accessory to HAAG-STREIT slit lamps and allows digital photography and video for documentation of the human eye. It is used at room temperature, usually by Ophthalmologists, Optometrists or Opticians in their consulting rooms, clinics, hospitals or teaching facilities.

Contents

1	Safety	5
2	Overview	
2.1	Camera Module CMX02 on BMV 900, BD 900	8
2.2	Camera Module CM 900 on BP 900	8
2.3	Release Module RMX01	8
2.4	Camera Module CMX02 on BQ 900	10
2.5	Release Module RMX01	10
3	Installation	
3.1	Install CMX02	13
3.2	Install cables	13
3.3	Image focus and position	13
3.4	Load compensation	15
3.5	Attaching the sticker for automatic left/right detection	15
3.6	Connecting the CMX02 and RMX01	15
4	Software / system starting up	17
5	System components	18
6	Accessories	19
7	Technical specifications	21

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1 Safety

Ambient conditions

•Transportation	Temperature	-40 °C	to	+70 °C
	Air pressure	500 hPa	to	1060 hPa
	Relative humidity	10%	to	95%
•Storage	Temperature	-10 °C	to	+55 °C
	Air pressure	700 hPa	to	1060 hPa
	Relative humidity	10%	to	95%
•Working	Temperature	+10 °C	to	+35 °C
	Air pressure	800 hPa	to	1060 hPa
	Relative humidity	30%	to	75%

Installation

- To avoid condensation, allow the instrument to adjust to room temperature for several hours before unpacking.
- Check that all components are fitted correctly (screwed connections, quick fit).
- All installations must be carried out by suitably trained personnel.

Operation and surrounding

- Only qualified and trained personnel should operate the equipment, the training is at the owner's responsibility.
- The image obtained via the eyepiece takes precedence when examining a patient.
- Use only HAAG-STREIT accessories.
- Shut down after every use. In case the dust cover is used: risk of overheating.
- Never operate the instrument in the same room with combustible gases, volatile solvents (alcohol, benzol) or flammable anesthetic agents.
- Avoid humidity.
- Check that the RMX01 is operating correctly before each examination.
- Make preliminary adjustments via the eyepiece rather than on-screen. This helps prevent the slit lamp from knocking against the patient's head.

Electrical

- Only hospital grade 3-conductor electrical power supply cables must be used. *For USA and Canada: Detachable Power Supply Cord Set, UL listed, type SJE, SJT or SJO, 3-conductor, not smaller than 18 AWG.*
- Plug, cable and ground lead connection of the socket have to be in perfect condition.
- The RM01 Release Module is secured with powerful magnets. Keep away from sensitive data-storage items.
- Use only PCs that conform to EN/IEC60601-1 standards, or operate them via a transformer that isolates them from the power supply. Generic EMC (electromagnetic compatibility) standards must be observed in accordance with EN 61000.
- PC peripherals (printers, monitors, scanners, etc.) must be operated via a transformer that isolates them from the power supply. These stipulations do not apply if a wireless connection (e.g. Bluetooth) is used.
- Ethernet connections must be configured with galvanic isolation in accordance with EN/IEC60601-1.
- Connect to the PC only with the USB-cables supplied.

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Environment

- Electric and electronic equipment must be separated from house waste! This equipment has been introduced into the market after the 13th of August 2005.
- Disposal via your local collecting point or your HAAG-STREIT distributor.
- Thus it is ensured, that no harmful materials get into the environment and that valuable raw materials can be used again.

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Cleaning

- Clean the housing only with a slightly water dampened cloth.
- No liquids, no alcohol or corrosive agents.
- Clean only the exterior surfaces of the optical parts.

Warranty / product liability

- The instrument should be operated in accordance with the chapter 'Safety'. Incorrect operation can damage the instrument. Thus no warranty claims can be accepted.
- Continued use of an instrument that has been damaged by incorrect operation can lead to personal injury. The manufacturer cannot accept liability in this case.
- Repairs and alterations on this equipment should only be carried out by HAAG-STREIT service technicians or by authorized persons.

Statutory requirements

- The Camera Module CM 900® has been designed to conform with the IEC / EN 60 601-1 and ISO 10939 standards. Manufacturing procedures, testing, commissioning, maintenance and repair are conducted under the observance of Swiss and international regulations.
- When combining different medical and/or nonmedical electrical equipment standard IEC / EN 60 601-1 applies.
- The Camera Module CM 900® must only be operated in situations where EN/IEC 60601-1 standards can be maintained, otherwise malfunctions may occur.
- The 'CE' marking confirms compliance of the Camera Module CM 900® with the directives 93/42/EEC and 89/336/EEC and the conformity module A.
- A copy of the declaration of conformity and information concerning EMC (Electromagnetic compatibility) of the present instrument can be requested at any time at HAAG-STREIT company.
- All statutory accident prevention regulations are to be observed.

Classification

IEC / EN Standard 60 601-1 Slit lamp accessory
safety class I
Application part Type B
Operation mode: continuous
operation

CE-Regulation 93/42 EEC Class I
FDA Class I

a


b

c

Pictograms

- a Read very carefully the instruction manual!
- b Disconnect the electrical supply plug!
- c Keep magnetic data-storage items away
RMX01 only:
- d Tip jack on RMX01 for power supply
- e Tip jack on RMX01 for CMX02
- f Tip jack on RMX01 for computer
- g Manufacturer
- h HS-Part Number
- i Type B equipment
- k Serial number
- l Date of manufacture

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
2 Overview

2.1 Camera Module CMX02 on BMV 900® or BD 900®

- 3 Fixation screw
- 4 Fixation ring
- 5 Camera Module CMX02 with logo (correct position of camera: opposite of microscope)
- 6 Background illumination cable
- 7 Cable holder
- 8 Camera cable
- 9 Connection cable computer (USB)
- 10 Release Module RMX01
- 11 Footswitch

HAAG-STREIT recommends the use of following components

- 1 Integrated background illumination (optical light lead)
- 2 Eyepiece with double crosshair reticule
- 14 Cold light source

2.2 Camera Module CMX02 on BP 900®

- 13 Fixation screw
- 15 Camera cable (USB)
- 16 Connection cable computer (USB)
- 22 Plug on camera head of Camera Module CMX02
- 23 Light diode (LED) indicator red = initialisation / error green = ready for use
- 24 Control lever beam splitter (70 / 30)
- 25 Fixation ring
- 27 Cable holder
- 28 Release Module RMX01

HAAG-STREIT recommends the use of following components

- 12 Background illumination pivoting mirror
- 26 Eyepiece with double crosshair reticule
- 14 Cold light source


2.3 Release Module RMX01

- 17 Right key = more exposure
- 18 Left key = less exposure
- 19 Trigger
- 20 Three stickers on lower surface
- 21 Remove protection foil and fix RMX01 on the base at right or left side of control lever

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2.4 Camera Module CMX02 on BQ 900

- 30 Cable holder
- 31 Camera cable (CMX02 - RMX01)
- 37 Light diode (LED) indicator
red = initialisation / error
green = ready for use
- 39 Control knob beam splitter
- 40 Fixation ring
- 41 Fixing screw at Video adapter


2.5 Release Module RMX01

- 32 Conduit tube
- 33 Cable tie
- 34 Computer cable (RMX01 - computer)
- 35 Power supply cable (RMX01 - power supply)
- 43 Sticker
left/right recognition
- 44 Key = more exposure
- 45 Key = less exposure
- 46 Trigger
- 47 Key = more exposure
- 48 Key = less exposure

- 36 Power supply with
- 50 Plug (country depending)

HAAG-STREIT recommends the use of following **components**

- 29 Background illumination
pivoting mirror
- 42 Eyepiece
with double crosshair reticule
- 49 Cold light source

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
3 Installation

3.1 Install CMX02

- Remove breath shield.
- Remove the cover on camera and adapter.
- Unscrew the 3 fixing screw (51) on the adapter. Tighten adapter (52) on camera.
- Rotate the camera into correct position and then affix it with 3 set screws (51). Setting image upright at PC.
- A distance screw (54) must be built in, thus the breath shield can be fixed correctly.

- Remove the cover.
- Align the positioning points (58) on the upper side of the parts to be joined.
- Push together and tighten firmly the joint in direction of the arrow by turning it.
- Position video adapter, tighten fixation ring (59).
- Loosen fixation screw (60) or (61) on video adapter. Tighten camera on adapter, rotate it into correct position, fix it with fixation screw.

3.2 Install cables

- Secure the camera cable (56) to the camera head (55), lead it eventually through the table and connect it to the computer. Possibilities: see chapter 2.
- On the slit lamp the cable is fixed on the microscope arm with two cable holder (57).

3.3 Image focus and position

Note the position of the camera! The camera must be mounted in the way, that the picture is correctly displayed on the screen. Rotate the camera in this way: the camera cable shows in direction of the microscope, so the HAAG-STREIT logo (53) is well visible.


In the case of misaligned images on the screen: loosen fixation screws (51) or (60) or (61), rotate the camera until the correct position image is displayed; fix camera again. The image can be inverted horizontally or vertically by the software.

Image focus on the monitor needs not to be adjusted, the interface has been designed so that optimal image focus is set automatically!

To achieve the same focus through both eye-pieces and monitor, the eye-pieces must be exactly adjusted with the focus test rod.

Or, the HAAG-STREIT recommendation: Correct focussing with eyepiece with double crosshair reticule.

When coupling the accommodation facility with convergent viewing, the normal refraction often has to be corrected by -1 to -2 dpt on the eye-piece.

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3.4 Load compensation (US Pat. 5.016.854)

The weight of components or accessories added to the microscope can be compensated by relieving springs so that height adjustment of the slit lamp remains smooth, no matter what accessories remain permanently mounted.

Compensating the weight of added components

- Turn the control lever down to its lowest stop, then loosen it by 1/4 turn.
- Swing the microscope and illumination unit to one side.
- Engage all three springs:

Engaging

Turn the screws (61) counter-clockwise until they are in the uppermost position.

Disengaging

Tighten the screws by turning them clockwise.

- Lift the microscope arm with one hand. Release it and check whether it returns to the original position. This will not happen unless the maximum load has been added.
- As a rule, it will be necessary to disengage as many relieving springs as required for ensuring that the arm and spring housing settle down in their original position.
- The load is adequately counter-balanced when the weight of illumination unit and microscope with all permanently mounted components exceeds the compensating force of the springs by a small amount.

3.5 Attaching the sticker for automatic left/right detection




- Detach the rail cover (66) and move the slit lamp out of the way. Clean the surface of the table.
- Remove the protective foil (62) of the sticker. Start with the corner opposite of the black surface.
- Position the sticker to the right-hand roller rail (64) and the edge of the table (65). Press down the black and white surface, eliminate any air bubbles.
- Carefully remove the rest of the 'Positionierungswerkzeug / Positioning tool' sticker (63) by tearing off along the perforation.
- Refit the slit lamp and rail cover.

3.6 Connecting the CMX02 and RMX01

IMPORTANT
ensure that all mains plugs are disconnected.

- Thread the computer cable (70) and power supply cable (68) into the conduit tube (67).
- Insert the angled plug (69) ...

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... of the power supply cable through the opening (74) at the back of the RMX01 release module, then plug it into the socket (77) on the opposite side.

- Press the cable into its cable holder (78), starting at the plug end. Press the cable into the groove (80) in the area near the control buttons.

- Push the computer cable (72) through the opening (74) in the RMX01, and plug it into the socket (77) marked with the PC symbol.

- Push the conduit tube with the computer and power supply cable through the opening (74) in the RMX01, and use a cable tie to secure it to the inside of the RMX01.

- Pull the conduit tube tight, and secure it at the other end with a second cable tie (84).

- Push the camera lead through the opening (74) in the RMX01, and plug it into the socket (76) marked with the camera symbol.

- Insert the RMX01 release module above the instrument base of the slit lamp. It is secured by four magnets (79). Model BQ 900® with screw (instrument base before 1998): remove first the 2 screws (83).

IMPORTANT

Keep magnetic data-storage items (such as credit cards) away from the magnets.

- Remove the protective foil from the cable holder (81) before sticking them to the microscope arm. Secure the camera lead (82).

- Connect the power supply (87) and plug (86), connect the computer cable (85) to the PC.

4 Software / system starting up

- Install the software supplied, as shown in the instruction manual.

After successful installation the **green** LED (88) will light up to indicate that the driver has been loaded, and




- the camera is ready for operation.

If the LED indicator lights up **red**

- check cables to ensure that they are connected


If the LED indicator does **not** light up

- check cables to ensure that they are connected
- switch on computer
- check power supply

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














For further information ask your HAAG-STREIT distributor

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7 Technical specifications

Minimum requirements for PC

Processor	Pentium IV
Clock frequency	min. 2 GHz for desktop min. 1.8 GHz for laptop
RAM	min. 1024 MByte, ideal 2048 MByte
Graphic card	PCI VGA card with DirectDraw support (not absolutely necessary, but you obtain best results)
Screen	Flat screen with min. 1280 x 1024 resolution, size ideally 19"
Interface	USB-2.0
Operating system	XP with SP2
Sound	best way is to work with sound

Key combinations

The following key combinations¹⁾ are used by the RMX01 / RMX01 (without F6 and F7) and may not be employed for other functions:

F6	Right eye
F7	Left eye
F10 / F11	Shorter exposure time
F9 / F12	Longer exposure time
F8	Trigger

¹⁾ in combination with EyeCap™ and EyeSuite™

Camera

Camera beam	model:	beam path:
	BD 900 BMV 900	left
	BP 900	HS recommends left
	BQ 900	HS recommends up, left or right

Sensor type	CCD
Interface	USB-2
Resolution	1600 x 1200 Pixel
Sensor dimension	1 / 1.8" (7.04 mm x 5.28 mm)
Framerate	12 fps (frames per second)

Object image section see Table

Circle: Object image section viewed through the microscope ocular.

Rectangle: Area of photograph. Indications on image section refer to dimensions of chip of 7.04 mm x 5.28 mm.

HAAG-STREIT Empfehlung (Gesichtsfeld) HAAG-STREIT recommandation (champ de vue) HAAG-STREIT recommendation (field of view)			
	Objektfeldgrösse ø mm Champ de l'objet ø mm Object field ø mm	Objektgrösse Kamera mm Dimension de l'objet caméra mm Object dimension camera mm	
		f 75	f 56
	12.5x Okular	4.76 x 5.95	4.76 x 5.95
	6.3x ø 32	16.7 x 20.8	22.3 x 27.9
	10x ø 20	10.5 x 13.2	14.1 x 17.6
	16x ø 12.7	6.7 x 8.3	8.9 x 11.2
	25x ø 8	4.2 x 5.2	5.6 x 7.0
	40x ø 5.1	2.7 x 3.3	3.6 x 4.5
	10x Okular	5.28 x 7.04	
	1x ø 18	17.9 x 23.9	BMV 900
	1.6x ø 11.3	11.2 x 14.9	BD 900

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Operation mode

continuous

Certified Quality Management System

EN ISO 13485
Reg. Nr. 11956

1250

HAAG-STREIT AG

Gartenstadtstrasse 10


3098 Koeniz, Switzerland


Phone ++ 41 31 978 01 11

Fax ++ 41 31 978 02 82

eMail info@haag-streit.com

internet www.haag-streit.com

HAAG-STREIT AG
Gartenstadtstrasse 10
3098 Koeniz, Switzerland
Phone ++ 41 31 978 01 11
Fax ++ 41 31 978 02 82
eMail info@haag-streit.com
internet www.haag-streit.com

11.09 - .05

1500.7220044.04020